



# Specification

## Acetic Acid 25%

### Reagent

#### Property

Assay  
Color  
Dilution Test  
Acetic Anhydride  
Chloride  
Residue After Evaporation  
Sulfate  
Titrable Base  
Reducing Substances Dichromate  
Reducing Substances Permanganate  
Heavy Metals  
Iron

#### Specification

24.0-26.0%  
10 APHA maximum  
To Pass Test  
0.01% maximum  
1 ppm maximum  
10 ppm maximum  
1 ppm maximum  
0.0004 meq/g maximum  
To Pass Test  
To Pass Test  
0.5 ppm maximum  
0.2 ppm maximum