



## Specification

### Acetone ACS Grade

<b><u>Property</u></b>	<b><u>Specification</u></b>
Assay	99.5% minimum
Color	10 APHA maximum
Aldehyde	0.002% maximum
IPA	0.05% maximum
Methanol	0.05% maximum
Residue After Evaporation	10 ppm maximum
Solubility	To Pass Test
Water	0.5% maximum
Titration Base	0.0006 meq/g maximum
Titration Acid	0.0003 meq/g maximum
Reducing Substances Permanganate	To Pass Test

Tested to: American Chemical Society Specifications (Tenth Edition) January 1, 2006