



## Specification

### Acetic Acid 50% Solution

#### PharmaPUR™

#### GMP Manufactured Product

**Property**

Appearance

Assay (w/w)

Sulfate

Chloride

Non Volatile Residue

Heavy Metals

Iron

Identification

**Specification**

Clear, colorless liquid; free from suspended matter

49.5-50.5%

Passes Test

25 ppm maximum

1 mg/20ml maximum

5 ppm maximum

5 ppm maximum

Passes Test

Made from Purified Water, USP and Glacial Acetic Acid, USP; which also meets EP and JP specifications.