



## Specification

### Nitric Acid 90% Fuming

### ChemPUR™ ACS Reagent

<u>Property</u>	<u>Specification</u>
Dissolved Oxides	0.1% maximum
Assay (w/w)	90% minimum
Sulfate	5.0 ppm maximum
Dilution Test	Passes Test
Residue after Ignition	20 ppm maximum
Chloride	0.7 ppm maximum
Arsenic	0.3 ppm maximum
Heavy Metals	5.0 ppm maximum
Iron	2 ppm maximum

Tested to: American Chemical Society Specifications (Tenth Edition) January 1, 2006