



2290 Avenue A, Bethlehem, PA 18017
E-mail: info@puritanproducts.com

PURITAN PRODUCTS
Custom Formulated and Specialty Chemicals

Phone 610-866-4225
FAX 610-866-4270

Specifications

SULFURIC ACID

20% Fuming

A.C.S. Reagent

<u>Property</u>	<u>Specification</u>
Appearance	Colorless to Very Light Brown
Assay (Free Sulfur Trioxide)	18.0 - 24.0%
MAXIMUM ALLOWABLE	
Residue after Ignition	0.002%
Nitrate (NO ₃)	1.0 ppm
Ammonium (NH ₄)	3.0 ppm
Arsenic (As)	0.03 ppm
Iron (Fe)	2.0 ppm

American Chemical Society Specifications, Reagent Chemicals, Tenth Edition, January 1, 2006