



Specification

Buffer Standard Solution pH 10.00 Blue

PharmaPUR™ NIST TRACEABLE

Property

Appearance
pH at 25 deg C

Specification

Clear, Blue Solution
10.00 +/- 0.01

Store at room or refrigerated temperatures

Our results are reported as described in Guide to Expression of Uncertainty in Measurement, ISO, Geneva, 1993, with a coverage factor of 2 (confidence level of 95%). The expanded uncertainty is $\pm 0.02\text{pH}$ for products in the pH range 0.50 through 10.00. Above 10.00pH, the expanded uncertainty is $\pm 0.05\text{pH}$.