



Specification

Hydrochloric Acid 3.00N Solution

PharmaPUR™ GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Assay Normality	2.95-3.05 N
Appearance	Clear, colorless liquid
Sulfate	Passes NF test
Residue after Ignition	2 mg maximum
Free Bromine or Chlorine	Passes NF test
Sulfite	Passes NF test
Heavy Metals	5 ppm maximum
Identification	Passes NF test

Made from USP Purified Water and Hydrochloric Acid, NF