



Specification

Hydrochloric Acid Solution 35.2 - 38.0%

ChemPUR™

Meets FCC Specifications

<u>Property</u>	<u>Specification</u>
Degrees Baume	To Pass Test
Appearance	Clear, colorless liquid
Assay (w/w)	35.2-38.0%
Specific Gravity	To Pass Test
Sulfate	5000 ppm maximum
Color (FCC)	Passes FCC Test
Extractable Organics*	5 ppm maximum
Oxidizing Substances	0.003% maximum
Non Volatile Residue	5000 ppm maximum
Reducing Substances	0.007% maximum
Iron	5 ppm maximum
Lead	1 ppm maximum
Mercury	100 ppb maximum
Identification	To Pass Test

Tested to: The current edition of Food Chemical Codex (Not a GMP manufactured product)