



## Specification

### Hydrogen Peroxide 29.0-32.0%

**ChemPUR™**

**ACS Reagent**

<u>Property</u>	<u>Specification</u>
Color	10 APHA maximum
Assay (w/w)	29.0-32.0%
Sulfate	5 ppm maximum
Ammonium	5 ppm maximum
Chloride	2 ppm maximum
Nitrate	2 ppm maximum
Phosphate	2 ppm maximum
Residue After Evaporation	20 ppm maximum
Titration Acid	0.0006 meq/g maximum
Heavy Metals	1 ppm maximum
Iron	0.5 ppm maximum

Tested to: American Chemical Society Specifications (Eleventh Edition) June 1, 2016