



Specification

Isopropyl Alcohol (2-Propanol)

ChemPUR™

ACS Reagent

<u>Property</u>	<u>Specification</u>
Color	10 APHA maximum
Assay	99.5% minimum
Carbonyl Compounds	0.002% maximum
Residue After Evaporation	10 ppm maximum
Solubility	To Pass Test
Water	0.2% maximum
Titration Acid or Base	0.0001 meq/g maximum

Tested to: American Chemical Society Specifications (Eleventh Edition) June 1, 2016