



Specification

n-Propyl Alcohol

ChemPUR™

Reagent

<u>Property</u>	<u>Specification</u>
Color	10 APHA maximum
Assay	99.5% minimum
IPA	0.05% maximum
Methanol	0.01% maximum
Ethyl Alcohol	0.01% maximum
Residue After Evaporation	10 ppm maximum
Solubility	Passes Test
Water	0.2% maximum
Titration Acid	0.0004 meq/g
Infrared Spectroscopy Scan	In accordance with the reference spectrum

Tested to: American Chemical Society Specifications (Eleventh Edition) June 1, 2016