



---

**Specification**  
**Sodium Hydroxide 10.0 N Solution**  
**PharmaPUR™**  
**GMP Manufactured Product**

<b><u>Property</u></b>	<b><u>Specification</u></b>
Assay Normality	9.90-10.10 N
Appearance	To Pass Test
Carbonate	0.6% maximum
Identification: Sodium <191>	Passes Test
Insoluble Substances and Organic Matter	To Pass Test
Chloride	30 ppm maximum
Heavy Metals	1 ppm maximum
Iron	2 ppm maximum
Endotoxin	2.0 EU/mL maximum

Made from USP Purified Water and Sodium Hydroxide, NF pellets which also meet EP and JP specifications.

The final solution is then filtered using a 0.2µm filter. The solution meets all requirements for residual solvents USP <467>.