



Specification

Sodium Hydroxide 1.00 N Solution

BioPUR™

<u>Property</u>	<u>Specification</u>
Assay Normality	0.95-1.05 N
Appearance	Clear and free from suspended matter
Color	10 APHA maximum
Carbonate	0.1% maximum
Sulfate	25 ppm maximum
Chloride	25 ppm maximum
Phosphate	5 ppm maximum
Nitrogen Compounds	0.001% maximum
Heavy Metals	15 ppm maximum
Iron	5 ppm maximum
Potassium	250 ppm maximum

Tested to: American Chemical Society Specifications (Tenth Edition) January 1, 2006