



Specification

Sodium Hydroxide 25% Solution

PharmaPUR™

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Appearance	To Pass Test
Carbonate	0.6% maximum
Assay (w/w)	24.50-25.50%
Identification: Sodium <191>	Passes Test
Insoluble Substances and Organic Matter	To Pass Test
Chloride	30 ppm maximum
Heavy Metals	1 ppm maximum
Iron	1 ppm maximum

Made from USP Purified Water and Sodium Hydroxide, NF pellets which also meet EP and JP specifications.

The final solution is then filtered using a 0.2µm filter. The solution meets all requirements for residual solvents USP <467>.