



Specification

Sodium Hydroxide 2.00 N Solution

ChemPUR™ Reagent

<u>Property</u>	<u>Specification</u>
Assay Normality	1.95-2.05 N
Appearance	Clear and free from suspended matter
Color	10 APHA maximum
Carbonate	0.1 % maximum
Sulfate	25 ppm maximum
Chloride	25 ppm maximum
Phosphate	5 ppm maximum
Nitrogen Compounds	0.001% maximum
Heavy Metals	15 ppm maximum
Iron	5 ppm maximum
Potassium	250 ppm maximum