



Specification

Sodium Hydroxide 5.0 N Solution

PharmaPUR™

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Assay Normality	4.9-5.1 N
Appearance	To Pass Test
Carbonate	0.6% maximum
Identification: Sodium <191>	Passes Test
Insoluble Substances and Organic Matter	To Pass Test
Chloride	30 ppm maximum
Heavy Metals	1 ppm maximum
Iron	0.5 ppm maximum
Endotoxin	2.0 EU/mL maximum

Made from USP Purified Water and Sodium Hydroxide, NF pellets which also meet EP and JP specifications.

The final solution is then filtered using a 0.2µm filter. The solution meets all requirements for residual solvents USP <467>.