



## Specification

### Sodium Hydroxide 50% Solution PharmaPUR™ GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Carbonate	3.0% maximum
Assay (w/w)	50.0-52.0%
Identification: Sodium <191>	Passes Test
Insoluble Substances and Organic Matter	Passes Test
Chloride	30 ppm maximum
Heavy Metals	1 ppm maximum
Iron	2 ppm maximum

Made from USP Purified Water and Sodium Hydroxide, NF pellets which also meet EP and JP specifications. The final solution is then filtered using a 2.5µm filter. The solution meets all requirements for residual solvents USP <467>.