



## Specification

### Nitric Acid 30%

#### Reagent Grade

<u>Property</u>	<u>Specification</u>
Appearance	Clear, colorless solution
Color	10 APHA maximum
Assay (w/w)	29.0-31.0%
Sulfate	1.0 ppm maximum
Residue after Ignition	5.0 ppm maximum
Chloride	0.5 ppm maximum
Arsenic	0.01 ppm maximum
Heavy Metals	0.2 ppm maximum
Iron	0.2 ppm maximum

Tested to: American Chemical Society Specifications (Eleventh Edition) June 1, 2016