



Specification

Sulfuric Acid 10.0 N Solution

PharmaPUR™

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Appearance of Solution	Clear and Colorless
Assay Normality	9.85-10.15 N
Identification A	To Pass Test
Identification B	To Pass Test
Limit of Nitrate	To Pass Test
Residue after Ignition	20 ppm maximum
Chloride	20 ppm maximum
Reducing Substances	To Pass Test
Arsenic	0.4 ppm maximum
Heavy Metals	2 ppm maximum
Iron	10 ppm maximum

Made from USP Purified water and Sulfuric Acid, NF. This product meets the USP guideline on residual solvents.