



Specification

Sulfuric Acid 1.00 N Solution

PharmaPUR™

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Assay Normality	0.95 - 1.05N
Appearance	Clear, colorless liquid
Residue after Ignition	50ppm maximum
Chloride	50ppm maximum
Reducing Substances	Passes NF test
Arsenic	1ppm maximum
Heavy Metals	5ppm maximum
Identification	Passes NF test

Made from USP Purified water and Sulfuric Acid, NF. This product meets the USP guideline on residual solvents.