



## Specification

### Sulfuric Acid 2.0 N Solution

#### PharmaPUR™

#### GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Description	Clear, colorless liquid
Assay Normality	1.95 - 2.05 N
Residue	0.005% maximum
Chloride	0.005% maximum
Reducing Substances	To Pass Test
Arsenic	1 ppm maximum
Heavy Metals	5 ppm maximum
Identification	To Pass Test

All raw materials meet the requirements stated by the USP/NF. Any water used to dilute said raw materials meets the USP requirements for Purified Water