



## Specification

### Sulfuric Acid 5.0 N Solution

#### PharmaPUR™

#### GMP Manufactured Product

| <u>Property</u>                    | <u>Specification</u> |
|------------------------------------|----------------------|
| Appearance of Solution             | Clear and Colorless  |
| Assay Normality                    | 4.9 - 5.1 N          |
| Identification A                   | To Pass Test         |
| Identification B                   | To Pass Test         |
| Limit of Nitrate                   | To Pass Test         |
| Residue after Ignition<br>Chloride | 20 ppm maximum       |
| Reducing Substances                | To Pass Test         |
| Arsenic                            | 0.4 ppm maximum      |
| Heavy Metals                       | 2 ppm maximum        |
| Iron                               | 10 ppm maximum       |

Made from USP Purified water and Sulfuric Acid, NF. This product meets the USP guideline on residual solvents.