



Specification

Conductivity Standard Solution 10 $\mu\text{S}/\text{cm}$

PharmaPUR™ NIST TRACEABLE

Property

Appearance
Conductivity @ 5 deg C
Conductivity @ 10 deg C
Conductivity @ 15 deg C
Conductivity @ 20 deg C
Conductivity @ 25 deg C
Conductivity @ 30 deg C
Conductivity @ 35 deg C

Specification

Clear, Colorless Solution
Report
Report
Report
Report
9.0-11.0 $\mu\text{S}/\text{cm}$
Report
Report

Store at room or refrigerated temperatures

Our results are reported as described in Guide to Expression of
Uncertainty in Measurement, ISO, Geneva, 1993, with a coverage factor of 2 (confidence level of 95%).