



Specification

Acetic Acid 50% w/v Solution

PharmaPUR™

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Appearance	Clear, colorless liquid; free from suspended matter
Assay (w/v %)	49.5 - 50.5%
Sulfate	Passes Test
Chloride	25 ppm maximum
Non Volatile Residue	1 mg/20ml maximum
Heavy Metals	5 ppm maximum
Iron	5 ppm maximum
Identification	Passes Test

Made from Purified Water, USP and Glacial Acetic Acid, USP; which also meets EP and JP specifications.