



Specification

Buffer pH 8.0

NIST TRACEABLE

<u>Property</u>	<u>Specification</u>
Appearance	Clear, Colorless Solution
pH at 5 deg C	Report
pH at 10 deg C	Report
pH at 15 deg C	Report
pH at 20 deg C	Report
pH at 25 deg C	8.00 +/- 0.08
pH at 30 deg C	Report

Store at room or refrigerated temperatures

Our results are reported as described in Guide to Expression of Uncertainty in Measurement, ISO, Geneva, 1993, with a coverage factor of 2 (confidence level of 95%). The expanded uncertainty is ± 0.02 pH for products in the pH range 0.50 through 10.00. Above 10.00pH, the expanded uncertainty is ± 0.05 pH.

Part Number 9743GMP