



Specification

Citric Acid 10%

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Description	Clear, colorless liquid
Assay	9.2-10.0% (w/w)
Assay (w/v %)	9.5-10.5%
Sulfate	0.015% maximum
Heavy Metals	0.001% maximum
Identification	Passes Test for Citrate

All raw materials meet the requirements stated by the USP/NF. Any water used to dilute said raw materials meets the USP requirements for Purified Water

Part Number 5049GMP