



Specification

Hydrochloric Acid 0.1N

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Description	Colorless liquid
Assay Normality	0.095-0.105 N
Identification	Passes NF test
Sulfate	Passes NF test
Residue after Ignition	2 mg maximum
Free Bromine or Chlorine	Passes NF test
Sulfite	Passes NF test
Heavy Metals	5 ppm maximum

Made from USP Purified Water and Hydrochloric Acid, NF

Part Number 6377GMP