



Specification

Sulfuric Acid 45%

GMP Manufactured Product

<u>Property</u>	<u>Specification</u>
Assay	44.0-46.0%
Identification	Passes Test
Residue after Ignition	50 ppm maximum
Substances Reducing Permanganate	Passes Test
Chloride	50 ppm maximum
Nitrate	10 ppm maximum
Arsenic	1 ppm maximum
Heavy Metals	5 ppm maximum
Iron	200 ppm maximum
Lead	5 ppm maximum

Part Number 2832GMP