



Specification

THF-Heptane Blend

<u>Property</u>	<u>Specification</u>
Appearance	Clear, colorless liquid
Color	10 APHA maximum
THF Assay	59-61%
Heptane Assay	39-41%
Residue After Evaporation	10 ppm maximum
Peroxidizing Substances	Not Detected by EM Quant Strips

Expiration Date: 12 months from peroxidizing compounds test date if drum remains unopened and is stored properly. Upon subsequent recertifications for peroxidizing compounds, under no circumstances will the cumulative expiration date of this product exceed 36 months

Part Number 8496